



Advanced Double Electric Breast Pump


Model Number 5161116

Designed with leading lactation professionals, the Evenflo Feeding Advanced Double Electric Breast Pump includes premium features moms are looking for. This hospital strength, full-time use pump helps maximize milk output, is comfortable and gentle on breasts, and is quiet and portable for discreet pumping anywhere.



BPA FREE

Parts intended to contact the breast or milk are made without BPA



HOSPITAL STRENGTH

- Effective, yet gentle
- Single and double pumping modes

MAXIMIZES MILK OUTPUT

- Independent speed and suction controls mimic baby's nursing pattern
- Multiple releases of milk

COMFORTABLE & GENTLE

- Includes 28.0mm and 30.5mm flanges
- Soft silicone insert

HYGIENIC & EASY TO CLEAN

- Closed system
- Prevents milk backup in the tubing

QUIET & PORTABLE FOR EVERYDAY USE

- 50% quieter and 30% lighter
- Three power sources: AC adapter, 6 AA batteries, car adapter*

FIRST CLASS CUSTOMER CARE

- Three year warranty
- One-on-one pumping support
- Access to breastfeeding education

*Sold Separately

Evenflo Feeding Advanced Double Electric Breast Pump

Model Number 5161116

WHAT'S INCLUDED

- Advanced Double Electric Breast Pump
- (2) 30.5mm Flange Kits
- (2) 28.0mm Soft Flange Inserts
- (2) 5 fl. oz. Advanced Milk Collection Bottles
- (2) Cap Rings & Sealing Discs
- (2) Extra Valves & Membranes
- (1) AC Adapter
- Balance + Standard Slow Flow Nipple*
- Instructions for Use
- Non-Retail Packaging



SPECIFICATIONS

ATTRIBUTE	ADVANCED DOUBLE ELECTRIC
Single/Double Pumping	Single/Double
Adjustable Suction Levels	Yes
Cycle Speed [Cycles per Minute]	30 - 80 (variable)
Overflow Protection	Yes (diaphragm)
Vacuum Range	Variable 25 - 250 mmHg
Cycling/Suction Control Mechanism	Microprocessor
Flange Nipple Tunnel Sizes [mm], Included with Pump	Medium (28.0mm) Soft Silicone Insert, Large (30.5mm) Hard Plastic Flange
Additional Flange Nipple Tunnel Sizes [mm], Sold Separately	X-Small (22.0mm) Soft Silicone Insert, Small (24.5mm) Soft Silicone Insert, X-Large (33.0mm) Soft Silicone Insert, XX-Large (36.0mm) Hard Plastic Flange
FDA Compliance	Class II: 510(k) K131153
Batteries	6 "AA" Batteries (not included)
Pump Unit Base	1.05 lbs.
Available Adapters, Sold Separately	Breast Milk Storage Bag Adapters Balance + Wide Pump Adapters

*The Balance + Slow Flow Nipple (0-6 months) is designed for breastfed babies and is compatible with the included Advanced Breast Milk Collection Bottles. It is also compatible with the Evenflo Feeding Balance + Standard Bottles (sold separately). The Balance + feeding system (bottles, nipples, pacifiers) was developed with Pediatric Feeding Specialists, and is recommended by 97% of parents.